

LAB NO. 755	HEATHCONTROL EQA SCHEME FOR DRUGS OF ABUSE IN URINE
EU workplace cut-offs	SAMPLE CODE NUMBER :- 103-1 February 2012

Your performance score = 17 out of 17 (100.0%)

QUANTITATIVE INTEGRITY TESTS	Group mean	SD	no.	Your data	SDs from mean	Score /out of	
creatinine (mmol/L)	5.40	0.31	56	5.57	0.55	1/1	
osmolality (mmol/kg)	342.00	3.80	10	No data			
pH	6.34	0.27	44	6.60	0.95		
specific gravity	1.010	0.002	15	No data			
TEST/DRUG/GROUP	Concentration present	Group report	Group %	no.	Your report	Errors	Score /out of
INTEGRITY TESTS							
aldehydes		negative	100.0	2			
nitrites		negative	100.0	2			
oxidants		negative	100.0	2			
sample check integrity		negative	100.0	10			
		OK	100.0	63	OK		1/1
DRUG/GROUP TESTS							
amfetamine group				71	present		
amfetamine	503 µg/L	present	97.7	43			
methyl-amfetamine	1502 µg/L	present	90.0	40			
cannabinoid group		not found	100.0	72	not found		1/1
benzoylecgonine		not found	100.0	70	not found		2/2
benzodiazepine group	660 µg/L	present	91.5	71	present		1/1
oxazepam	660 µg/L	present	100.0	24			
methadone		not found	100.0	57	not found		2/2
EDDP		not found	100.0	36	not found		2/2
buprenorphine+metab		not found	100.0	43	not found		1/1
opiate group	627 µg/L	present	93.0	71	present		2/2
6-monoacetylmorphine	21 µg/L	present	79.5	39	present		4/4
morphine	606 µg/L	present	97.7	43			
others				42			
barbiturate group		not found	100.0	21	not found		



EINGEGANGEN
26. April 2012
TQM

LAB NO. 755	HEATHCONTROL EQA SCHEME FOR DRUGS OF ABUSE IN URINE
EU workplace cut-offs	SAMPLE CODE NUMBER :- 103-2 February 2012

Your performance score = 14 out of 14 (100.0%)

QUANTITATIVE INTEGRITY TESTS				Group mean	SD	no.	Your data	SDs from mean	Score /out of
creatinine (mmol/L)				1.53	0.15	58	1.39	-0.91	1/1
osmolality (mmol/kg)				84.80	4.64	10	No data		
pH				7.58	0.48	43	7.50	-0.16	
specific gravity				1.003	0.002	16	No data		
TEST/DRUG/GROUP	Concentration present	report	Group %	no.	Your report	Errors	Score /out of		
INTEGRITY TESTS									
aldehydes		negative	100.0	2					
nitrites		negative	100.0	2					
oxidants		negative	100.0	2					
sample check		negative	100.0	9					
integrity	creatinine low	not OK	79.0	62	not OK		1/1		
DRUG/GROUP TESTS									
amfetamine group		not found	100.0	66	not found		1/1		
cannabinoid group		not found	100.0	67	not found		1/1		
benzoylecgonine		not found	100.0	66	not found		2/2		
benzodiazepine group				66	not found				
methadone		not found	100.0	54	not found		2/2		
EDDP		not found	100.0	32	not found		2/2		
buprenorphine+metab		not found	100.0	43	not found		1/1		
opiate group		not found	100.0	66	not found		1/1		
6-monoacetylmorphine		not found	100.0	25	not found		2/2		
others				41					
barbiturate group		not found	100.0	18	not found				

LAB NO. 755	HEATHCONTROL EQA SCHEME FOR DRUGS OF ABUSE IN URINE
EU workplace cut-offs	SAMPLE CODE NUMBER :- 103-3 February 2012

Your performance score = 15 out of 15 (100.0%)

QUANTITATIVE INTEGRITY TESTS				Group mean	SD	no.	Your data	SDs from mean	Score /out of
creatinine (mmol/L)				9.18	0.54	56	9.43	0.47	1/1
osmolality (mmol/kg)				564.89	5.56	9	No data		
pH				5.98	0.24	41	6.10	0.52	
specific gravity				1.017	0.003	14	No data		
TEST/DRUG/GROUP	Concentration present	Group report	Group %	no.	Your report	Errors	Score /out of		
INTEGRITY TESTS									
aldehydes		negative	100.0	2					
nitrites		negative	100.0	2					
oxidants		negative	100.0	2					
sample check integrity		negative	90.0	10					
		OK	98.3	59		OK	1/1		
DRUG/GROUP TESTS									
amfetamine group		not found	82.9	70	not found		1/1		
cannabinoid group		not found	100.0	71	not found		1/1		
benzoylecgonine		not found	100.0	70	not found		2/2		
benzodiazepine group		not found	100.0	70	not found		1/1		
methadone		not found	100.0	57	not found		2/2		
EDDP		not found	100.0	33	not found		2/2		
buprenorphine+metab		not found	100.0	43	not found		1/1		
opiate group	3751 µg/L	present	95.8	71	present		1/1		
6-monoacetylmorphine		not found	100.0	39	not found		2/2		
morphine	751 µg/L	present	48.8	43					
morphine glucuronide	751 µg/L	present	100.0	1					
codeine	2999 µg/L	present	100.0	43					
others				45					
barbiturate group		not found	100.0	21	not found				